




PATIENT HANDLING

ALUMINIUM SCOOP STRETCHER

Aluminium stretcher with break-apart scissor action at both ends



Ergonomic design allows for ideal immobilisation during patient transfer.

PRODUCT FEATURES

- Lightest scoop on the market weighing only 16 pounds
- Superior Comfort and safety to patients who need spinal immobilization
- Eliminates the need for log-roll maneuvers, decreasing movement to cervical spine.
- Two hinged, interlocking pieces allow operators to bring the two halves together beneath the patient and gently scoop them up.
- 159kg (350lbs) load capacity
- Easy to clean
- Adjustable length which can accommodate multiple- height patients
- Folding design which provides compact storage in ambulance compartments.

SPECIFICATION (OPEN SIZE)

Width	55cm (21.6 inches)
Depth	12cm (4.7 inches)
Length	214cm (84.2 inches)
Weight	7.26kg (16 lbs)
Weight Capacity	159kg (350 lbs)

SPECIFICATION (FOLDED SIZE)

Width	20.9cm (8.2 inches)
Depth	14cm (5.5 inches)
Length	59cm (22.8 inches)
Weight	7.26kg (16 lbs)
Weight Capacity	N/A