



### PATIENT HANDLING

## BARIATRIC MATTRESS

Designed for evacuating obese people up to 300kg (661lbs) in lying position



The Bariatric Mattress needs to be operated by several trained persons. Also suitable for spiral staircases.

The Bariatric Mattress can be used to move evacuees along the floor (by using at least 4 persons and using the additional towing straps) and have 12 handles placed on the sides to allow lifting and easy maneuverability.

### KEY BENEFITS

- TÜV certified on 300kg (661lbs)
- Extra wide version
- Made of high quality materials (fire resistance classification M2)
- Extra-long towing straps to help operators with improved handling practices
- Handgrips to assist with sideways manoeuvrability
- Carrying people is no longer necessary
- Ready for use within seconds.
- Safe evacuation down the stairs.
- Speed of descent is determined by the user.
- Fits into your emergency evacuation plan
- CE, TÜV NORD and ISO tested.

### SPECIFICATION

|  |                     |
|--|---------------------|
| Length   | 200cm (78.7 inches) |
| Width  | 86cm (33.8 inches)  |
| Thickness  | 1.3cm (0.5 inches)  |
| Weight   | 10kg (22.1 lbs)     |
| Maximum payload                                    | 300kg (661 lbs)     |
| Speed of descent to be determined by the user      | ✓                   |
| Heat and flame retardant                           | ✓                   |
| Fire resistance classification M2                  | ✓                   |
| Blood, urine, acid & alkaline resistant            | ✓                   |
| Equipped with headrest, body covers and foot cover | ✓                   |
| Ready to use in seconds                            | ✓                   |
| Manufacturers warranty                             | ✓                   |