



### PATIENT HANDLING

## MULTI-PURPOSE EMS COT MATTRESS

The complete solution for trauma patient transportation



The multi-purpose **EMS Cot Mattress** is the perfect solution for multifaceted rescue needs that has multiple features in 1:

- A stretcher function
- A dragging function
- A carry function

The multi-purpose **EMS Cot Mattress** has both towing straps and side handles so it can be used directly as an evacuation AND rescue mattress. The great advantage is that the mattress fits straight back on the stretcher, ready for transport. The multi-purpose **EMS Cot Mattress** has been designed to fit comfortably on Stryker, Ferno and other stretchers.

### KEY BENEFITS

- Available to fit any size stretcher
- Adjustable foot cover (patented)
- Made of high quality & fire resistant materials
- Blood, urine, acid and alkali resistant
- OHS & ergonomic accountable.
- No more unnecessary transfers needed
- Ready to use within seconds
- Cost and space saving in the ambulance.
- CE & TÜV certified.

### SPECIFICATION

Length	180.4cm/71"
Width	47cm/18.5"
Thickness	8.1cm/3.2"
Weight	11kg/24.3 lbs
Multi-purpose solution for the professional user	✓
Ergonomic friendly	✓
Available for any size stretcher	✓
Fire resistance classification M2	✓
Blood, urine, acid & alkaline resistant	✓
Equipped with headrest, body covers and foot cover	✓
Ready to use in seconds	✓
Manufacturers warranty	✓

#### Breakaway sections



#### Reinforced bottom



#### Head support



#### Adjustable foot cover

