



PATIENT RESCUE

### LIFEPAD

Provide first aid safely and quickly –  
Your helper in an emergency



Most cardiac arrests happen at home. Many first-aiders are initially overwhelmed when confronted with an emergency. The LifePad® guides you through the resuscitation process and helps you to perform the right steps properly.

- Visual feedback through coloured LEDs to maintain the optimal compression pressure
- Acoustic signal to specify the optimal compression rate of 100 beats/minute
- Easy positioning thanks to instructions and marking directly on the pad

#### PRODUCT FEATURES

- Portable and easy to store thanks to compact size and flexible material
- Raised surface part for the hands helps with correct positioning
- Adapts to the individual's body shape
- Can be used in combination with a defibrillator (sound can be switched off)
- Patented technology
- Can be used by and on persons aged 12 and over
- IP 44 splash-proof
- Medical device
- Product dimensions (LxWxH): 210 x 160 x 3 mm (Height Housing: 15 mm)
- Product weight: approx. 130 g
- 5-year warranty

