



PEDIATRIC TRANSPORT

PEDIATRIC VACUUM STRETCHER

Pediatric vacuum mattress with integrated harness



The **Paediatric Vacuum Mattress** with integrated sling is intended for the safe transport of newborns and children. It is therefore the most complete device for adapted paediatric transport and immobilisation.

The **Pediatric Vacuum Mattress** is quickly attached to the stretcher with 5 black attachment straps. Each strap is adjustable and has quick release buckles. The red base of the **Pediatric Vacuum Mattress** is quickly attached to the stretcher with 5 attachment straps. Each strap is adjustable and has a metal safety buckle, which is fitted directly to the stretcher frame.

To facilitate installation, each type of strap has its own colour. This makes it easy to identify the type of strap used on the stretcher: black for the straps used to attach the device to the stretcher and blue for the holding straps.

FEATURES

- For the safe transport of a child on an ambulance stretcher.
- From newborns to children aged 6 to 7 years (125m).
- Quick attachment to the stretcher.
- Optional carry bag.
- In red/blue polyurethane.
- Dimensions : 134.6cm x 60cm (53 x 23.6 inches).
- Weight : 3kg (6.6lbs).
- 5 coloured polypropylene straps.
- 2 support straps and 2 harness straps.
- Guarantee : 5 years.

SPECIFICATION

Width	60cm (23.6 inches)
Depth	4.5cm (1.77 inches)
Length	134.6 (53 inches)
Weight	3kg (6.6 lbs)