



### PATIENT HANDLING

## PRO-EMS CHAIR HYBRID

Combines the benefits of a carry chair and also allows for an attachable powered track system



Combines the benefits of a carry chair and also allows for an attachable powered track system which enables patient transportation up and down stairs. The attachable track system minimises the weight of the chair so it can easily be carried by operators.

### PRODUCT FEATURES

- Super lightweight hybrid power
- Can be used as a carry chair or power chair
- 4-speed control panel
- Simple to attach track power system
- Height adjustable rear handle
- Footrest
- Easy to store in EMS vehicles

### SPECIFICATION

Height	95.5 cm (37.6 inches)
Width	52 cm (20.5 inches)
Depth	20.3 cm (8 inches)
Weight	25.6 kg (56.50lbs)
Weight Capacity	254 kg (560lbs)

