



PATIENT HANDLING

PRO-EMS CHAIR STANDARD

The lightest EMS chair with the highest weight capacity



The **PRO-EMS-S** is the market's lightest EMS Chair with the highest weight capacity. The **PRO-EMS-S** simplifies patient transportation with its single person operation and track-belts which give EMS professionals ultimate control for descending transportation.

PRODUCT FEATURES

- Super lightweight
- Anti-slip height adjustable handle
- Additional operator handles
- Quick release 4-point harness
- Easy clean ABS seating
- Industrial locking rear wheels
- Passenger footrest

SPECIFICATION

Height	95.5 cm (37.6 inches)
Width	52 cm (20.5 inches)
Depth	20.3 cm (8 inches)
Weight	14.2 kg (31.3lbs)
Weight Capacity	254 kg (560lbs)

