



PATIENT RESCUE

SCOOP STRETCHER

Plastic stretcher with break-apart scissor action at both ends



Ergonomic design allows for ideal immobilisation during patient transfer.

PRODUCT FEATURES

- Lightest scoop on the market weighing only 7.3kg (16lbs) 1kg (2lbs) less than major brand name)
- Maximum Length is 234cm (92") a full 33cm (13") longer than competitors
- Superior Comfort and safety to patients who need spinal immobilization
- Eliminates the need for log-roll maneuvers, decreasing movement to cervical spine.
- Two hinged, interlocking pieces allow operators to bring the two halves together beneath the patient and gently scoop them up.
- Made from lightweight, X-ray translucent, high impact, composite materials and includes three restraints.
- 227kg (500lbs) load capacity
- Easy to clean
- Adjustable length which can accommodate multiple- height patients
- Folding design which provides compact storage in ambulance compartments.

SPECIFICATION

Width	90.2cm (35.5 inches)
Depth	10.2cm (4 inches)
Maximum Length	234cm (92 inches)
Minimum Length	165cm (65 inches)
Weight	7.25kg (16 lbs)
Weight Capacity	227kg (500 lbs)

