



PATIENT RESCUE

### VIVEST X-SERIES AED

Designed with mobility, robustness, ease of use and reliability in mind



Slimline and lightweight means that the ViVest X-Series AED can fit into most bags and be transported to various locations.

The X1 AED has clear visual images, enhanced further using LED lighting and loud verbal commands, while the X3 has on screen instructions to guide the user and provide the patient with the best chance of survival.

Daily self testing, drop tested to 1.5 metres, IP55 protected without the need for a carry case and a fail-safe button ensuring the AED is always ready to save lives.

#### PRODUCT HIGHLIGHTS

- Slimline and lightweight – H 23.2 x W 20.9 x D 5.9cm, 1.5kg
- One of the easiest to use on the market
- Time to shock 7 seconds
- Battery life of 5 years and capable of up to 200 shocks
- Smart Pads life of 5 years which allow for biphasic shock
- CPR features with metronome
- Algorithm capable of detecting implants such as pacemaker and internal defibrillator
- Automated ECG recording of up to 8 hours
- Daily self testing ensures device is ready for use
- Drop tested to 1.5 metres
- Shock button has built in fail-safe button to alleviate risk of failure
- Able to identify a pulse and make the decision not to shock

#### SPECIFICATION

Size (H x W x D)	23.2 x 20.9 x 5.9cm (9.16 x 8.22 x 2.3 inches)
Colour display	X3
Multi-language option	X3
Protection grade	IP55
Daily self check	Yes
Battery life	5 Years

