



PEDIATRIC TRANSPORT

EMERGENCY CHILD RESTRAINT

For the safe and effective transport of infants and children in an ambulance



SMALL

4-18lb
2-8kg

MEDIUM

13-46lb
6-21kg

LARGE

42-110lb
19-50kg

The revolutionary **Emergency Child Restraint** is designed to safely secure children during EMS transport in any emergency vehicle. Since it is not acceptable for a child to be transported unrestrained on a cot, whether it is a new-born or even a teenager, the **Emergency Child Restraint** provides the perfect solution.

The 3 color coded **Emergency Child Restraint** harnesses covers a weight range of between 4 lbs (1.8kg) and 110 lbs (50kg), offering the maximum protection available for children of any size and weight. **The Emergency Child Restraint** has been crash tested under the strictest of environments. Its design has no stitching or velcro providing the ultimate protection and safety for you and your patient in every use/application.

Designed with patient comfort in mind, the **Emergency Child Restraint** has a unique soft touch system, with full clinical access to the torso, without removing the restraint during transport.

FEATURES

- Transport children from 4lb (1.8 kg) to 110 lbs (50 kg).
- 3 color coded integrated, fully adjustable **Emergency Child Restraint** harnesses.
- Includes 4 extra-long universal cot straps for any cot.
- Ergonomically designed to fit children of all sizes/shapes.
- Unique quick reference color coded labelling system in line with industry emergency treatment of children/ color dosing systems.
- Full clinical access to a child's torso for medical intervention, without removing restraint while in use.
- Special fixing system, ensuring no downwards pressure on infant during its use.
- Uniquely designed and protected soft touch system, which reduces vibration for ultimate infant comfort.
- No stitching or velcro, ensuring maximum infection control.
- Easy to maintain bio safe cot straps and restraint material ensuring 100% cleanability.
- Fully Crash tested.

