



PATIENT RESCUE

### VACUUM MATTRESS

Used for the immobilization of patients, especially in the case of vertebra, pelvis or limb trauma



This flexible pad allows a perfect moulding of the patient, and can be folded into three. The mattress consists of 3 canvases, assembled by high frequency welding. The upper, middle and lower canvases are made of PVC .

#### PRODUCT FEATURES

- Reduces pressure and provides better support for lumbar injuries and spinal trauma
- Minimizes body movement compared to spinal backboards
- Protects against cold in inclement weather conditions
- Significantly superior stability and comfort compared to a spinal backboard
- Available in adult and pediatric sizes
- Packs away into a carry bag & supplied with pump
- Washable & X-ray translucent
- Fire retardant
- Waterproof
- Can be used in temperatures from -40o to +60o (140oF)

#### SPECIFICATION

Width	90.2cm (35.5 inches)
Depth	7cm (2.75 inches)
Length	221cm (86.7 inches)
Weight	8.2kg (18 lbs)
Weight Capacity	150kg (330 lbs)